

PREPARATION

I. SITE PREPARATION

This section of the manual describes how to prepare the site for installing your UNcle and how to get the instrument up and running. There are important safety warnings included here which must be observed by all personnel that have access to the instrument. Hardware installation should be completed before installing the software.

Do not unpack the instrument until Unchained Labs personnel is present. An Unchained Labs qualified engineer or field applications scientist should provide assistance with the site preparation and installation of your UNcle. The following section, therefore, is for information only.

Before installing your UNcle in the workplace, preparation of the site according to the guidelines in this section is required. At the end of the section there are a set of checklists to help ensure that all the preparation tasks have been addressed.



IMPORTANT - If your UNcle Instrument is moved after it has been installed, it must be re-qualified by an Unchained Labs service engineer before use.

System Components

Your UNcle consists of the following items:

- UNcle system
- External computer system with touch screen monitor, keyboard and mouse
- Connector cables

Dimensions and Weights

The dimensions and weights of a typical UNcle system are given below. Ensure that the installation site can accommodate the dimensions and weights.

Component	Width	Depth	Height	Weight
UNcle	54 cm (21.3 in)	50 cm (20 in)	58 cm (22.8 in)	50 kg (110 lb)
Desktop computer	17 cm* (7 in)	43 cm* (17 in)	37 cm* (15 in)	10 kg* (22 lb)

Clearances

Once in operation, your UNcle system needs at least 5 cm (2 inches) clearance on all sides for ventilation. Do not obstruct airflow to the air inlet grill on the bottom panel. For example, do not remove the feet so the instrument is resting directly on its base.

For installation and on-site repair and maintenance, the engineer will require at least 60 cm (24 inches) clearance either side of the instrument to gain access to the side panels.

Decontamination



WARNING – Biohazard. It is your responsibility to decontaminate components of the instrument before requesting service by a service engineer. Unchained Labs will not accept items that have not been decontaminated where it is applicable to do so.

Cleaning your UNcle

Always turn the power switch off and disconnect the power cord from the main power source before using liquids to clean the instrument.

- Wipe up all spills immediately.
- Periodically clean the outside surfaces of the instrument using a cloth or sponge that has been lightly dampened with water.
- If required, clean the surfaces using a mild soap solution diluted with water and then wipe with a damp cloth or sponge to remove all residue.
- If needed, clean the stage using a cloth or sponge that has been lightly dampened with water.
- If a bleach solution has been used, wipe the instrument using a lint-free cloth that has been lightly dampened with water to remove the bleach residue.

CAUTION! Do not use abrasive cleaners. Do not spray cleaner directly onto the instrument. Do not let water or other fluids drip inside the instrument.

CAUTION! If there are any questions or concerns about decontamination or cleaning, please contact your local Unchained Labs service provider.

Environmental Requirements

The site should be maintained under the following conditions.

The user should be aware that fluctuations in environmental conditions can cause unexpected changes in the instrument's performance. For best results, ensure that the temperature and humidity are maintained at a constant level.

Condition	Acceptable Range
Operating Temperature	18 to 28 °C (64 to 83°F)
Storage Temperature	0 to 45 °C (32 to 113°F)
Humidity	40 to 80 % relative humidity, non-condensing, altitude <2000m, sound level <80db
Pollution degree	2
Overvoltage category	II
Mains voltage fluctuation	+/- 10%

This apparatus is for **indoor use only**. Under normal operation it has ingress protection rating IP20. The instrument should not be placed in direct sunlight.

Avoid significant and sudden changes to the environmental conditions immediately preceding use of the instrument; for example when unpacking, installing or relocating your UNcle system.

Electrical Requirements

Your UNcle requires connection to one standard electrical outlet.

The desktop computer requires connection to one electrical outlet. If there is a separate monitor then this may require connection to an additional outlet (or it may plug into the desktop computer). The desktop computer will have its own standard power requirements.



IMPORTANT - In case of emergency, the main power supply to the instrument must be disconnected immediately. Ensure that the UNcle location allows for this.

Your UNcle is factory configured to operate interchangeably at ~120 VAC or ~240 VAC.

Voltage	Frequency	Fuse Rating	Max. Power
220 – 240 V	50 – 60 Hz	T6.3AL,250V anti-surge	600 W (VA)
100 – 125 V			

Over current protection >110%.

Checklists

Check that the location fulfils the following criteria:

- ☒ The bench or table area is sufficiently large and can withstand the weight of the instrument.
- ☒ There are sufficient electrical outlets within range.
- ☒ The location of the instrument and computer allows for appropriate ergonomic usage and access to the power connection for safety.

Check that all the following components are present:

- ☒ UNcle system
- ☒ Power lead x 3
- ☒ USB lead
- ☒ All-in-One Workstation
- ☒ Mouse
- ☒ Keyboard

If the installation site is suitable and all of the components are present, installation can commence. The installation engineer will set up the instrument, giving the user the opportunity to read the important safety information in the following section, so that potential hazards can be spotted when operation of the instrument begins.